



Product Datasheet

MafA Antibody (orb1529912)


Description	MafA Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 300 and C-term
Target	MafA
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8CF90
NCBI	NP_919331.1
Dilution Range	WB - 1:1,000 - 1:10,000
Expiration Date	12 months from date of receipt.



Detection of MafA, MafB and cMaf by west...



Samples: Immunoprecipitated cross-linked...



Electrophoretic Mobility Shift of MafA, ...