



Product Datasheet

BTF Antibody (orb1529917)


Description	BTF Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 150 and 200
Target	BTF
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NYF8
NCBI	NP_055554.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 1 - 10 µg/mg lysate; IHC - 1:500 to 1:2,000. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human BTF by western blot.S...



Detection of mouse BTF by western blot.S...



Detection of human BTF by immunohistoche...