


Product Datasheet

BID, Phospho (S78) Antibody (orb1530051)

Description	BID, Phospho (S78) Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	surrounding serine 78
Target	BID, Phospho (S78)
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	500 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P70444
NCBI	NP_031570.2
Dilution Range	WB - 1:250 - 1:1,000
Expiration Date	12 months from date of receipt.



Detection of mouse Phospho-BID (S78) aft...