


Product Datasheet

p70S6K Antibody (orb1530085)

Description	p70S6K Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 475 and C-term
Target	p70S6K
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P23443
NCBI	NP_003152.1
Dilution Range	WB - 1:2,500 - 1:20,000
Expiration Date	12 months from date of receipt.



Detection of human and mouse p70S6K by W...