


Product Datasheet

KPNB1 Antibody (orb1530131)

Description	KPNB1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 825 and C-term
Target	KPNB1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q14974
NCBI	NP_002256.2
Dilution Range	WB - 1:5,000 - 1:15,000; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of KPNB1 by WB and IP. Samples:...