


Product Datasheet

MnSOD Antibody (orb1530188)

Description	MnSOD Antibody
Species/Host	Goat
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 75 and 125
Target	MnSOD
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P04179
NCBI	NP_000627.2
Dilution Range	WB - 1:1,000 - 1:4,000
Expiration Date	12 months from date of receipt.



Detection of
MnSOD by
western
blot.Sampl...