



Product Datasheet

PNUTS Antibody (orb1530198)

Description	PNUTS Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 900 and C-term
Target	PNUTS
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q96QC0
NCBI	NP_002705.2
Dilution Range	WB - 1:1,000 - 1:5,000; IP - 1 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human PNUTS by western blot...



Detection of human PNUTS by WB and IP.Sa...