



Product Datasheet

DBC1/p30 DBC Antibody (orb1530206)


Description	DBC1/p30 DBC Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 875 and C-term
Target	DBC1/p30 DBC
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8N163
NCBI	Q6P0Q9
Dilution Range	WB - 1:1,000 - 1:10,000; IP - 1 - 10 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of mouse DBC1 by western blot....



Detection of human DBC1/p30 DBC by immun...



Detection of human DBC1/p30 DBC by WB an...