


Product Datasheet

CRSP3 Antibody (orb1530221)

Description	CRSP3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 1300 and C-term
Target	CRSP3
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9ULK4
NCBI	NP_004821.2
Dilution Range	WB - 1:1,000 - 1:5,000; IP - Not recommended. Rabbit anti-CRSP3 Antibodies A300-423A and A300-424A are recommended for IP.
Expiration Date	12 months from date of receipt.



Detection of human CRSP3 by western blot...