


Product Datasheet

CRSP3 Antibody (orb1530223)

Description	CRSP3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 400 and 450
Target	CRSP3
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9ULK4
NCBI	NP_004821.2
Dilution Range	WB - Not recommended. Use rabbit anti-CRSP3 antibody A300-425A.; IP - 1 - 4 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of
human CRSP3
by
immunoprecip...