



Product Datasheet

EWS Antibody (orb1530230)


Description	EWS Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 100 and 150
Target	EWS
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q01844
NCBI	NP_005234.1
Dilution Range	WB - 1:10,000 - 1:20,000; IP - 1 - 5 µg/mg lysate; IHC - 1:2,000 - 1:10,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human and mouse EWS by west...



Detection of human EWS by immunohistoche...



Detection of human EWS by western blot o...