


Product Datasheet

EWS Antibody (orb1530232)

Description	EWS Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 1 and 50
Target	EWS
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q01844
NCBI	NP_005234.1
Dilution Range	WB - 1:5,000 - 1:15,000; IP - Not recommended. Rabbit anti-EWS Antibodies A300-417A and A300-418A are recommended.
Expiration Date	12 months from date of receipt.



Detection of human EWS by western blot.S...