



Product Datasheet

DHX15 Antibody (orb1530270)

Description	DHX15 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	between 350 and 400
Target	DHX15
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O43143
NCBI	NP_001349.2
Dilution Range	WB - 1:5,000 - 1:15,000; IP - Not recommended. A300-390A has performed well and is suggested for IP.; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human DHX15 by western blot...



Detection of human DHX15 by immunohistoc...