



Product Datasheet

Bcl11b Antibody (orb1530279)


Description	Bcl11b Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 500 and 550
Target	Bcl11b
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9C0K0
NCBI	NP_612808.1
Dilution Range	WB - 1:5,000 - 1:15,000; IP - 1 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human Bcl11b by western blo...



Detection of human Bcl11b by immunohisto...



Detection of human Bcl11b by WB and IP.S...