



Product Datasheet

MLL1 Antibody (orb1530291)


Description	MLL1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, IP, WB
Immunogen	between 2725 and 2775
Target	MLL1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q03164
NCBI	NP_005924.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 -5 µg/mg lysate; ChIP - 1 - 5 µg. Previous lots of this antibody have performed in this application.
Expiration Date	12 months from date of receipt.



Detection of human MLL1 by western blot....



Detection of human MLL1 by western blot ...



ChIP-chip scatter plot of anti-MLL1 enrichment...