



Product Datasheet

NCOA2/SRC2 Antibody (orb1530327)


Description	NCOA2/SRC2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 1400 and C-term
Target	NCOA2/SRC2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q15596
NCBI	NP_006531.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 1 - 4 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of mouse NCOA2/SRC2 by wester...



Detection of human NCOA2/SRC2 by immuno...



Detection of human NCOA2/SRC2 by WB and ...