



Product Datasheet

CoAA Antibody (orb1530354)


Description	CoAA Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	between 150 and 200
Target	CoAA
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q96PK6
NCBI	NP_006319.1
Dilution Range	WB - 1:1,000 - 1:10,000; IHC - 1:250 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human CoAA by immunohistochemistry.



Detection of human CoAA by western blot.



Detection of human CoAA by western blot.