



Product Datasheet

PHD4 Antibody (orb1530356)

Description	PHD4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 450 and C-term
Target	PHD4
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NXG6
NCBI	NP_808808.1
Dilution Range	WB - 1:500 - 1:2,500
Expiration Date	12 months from date of receipt.



Detection of human PHD4 by western blot....



Detection of human PHD4 by western blot....