


Product Datasheet

FUS Antibody (orb1530403)

Description	FUS Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 200 and 250
Target	FUS
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P35637
NCBI	NP_004951.1
Dilution Range	WB - 1:200 - 1:10,000; IP - Not recommended. Rabbit anti-FUS Antibodies A300-302A, A300-293A and A300-294A are much more efficient in IP.
Expiration Date	12 months from date of receipt.



Detection of human and mouse FUS by west...