



Product Datasheet

Claspin Antibody (orb1530449)


Description	Claspin Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 450 and 500
Target	Claspin
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9HAW4
NCBI	AAG24515.1
Dilution Range	WB - 1:500 - 1:2,500; IP - 1 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human Claspin by western bl...



Detection of human Claspin by immunohist...



Detection of human Claspin by western bl...