


Product Datasheet

Claspin Antibody (orb1530452)

Description	Claspin Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 450 and 500
Target	Claspin
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9HAW4
NCBI	AAG24515.1
Dilution Range	WB - 1:1,000 - 1:5,000; IP - Not recommended. Use rabbit anti-Claspin antibody A300-266A or A300-267A.
Expiration Date	12 months from date of receipt.



Detection of human Claspin by western bl...