


Product Datasheet

RPA32, Phospho (S33) Antibody (orb1530476)

Description	RPA32, Phospho (S33) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	surrounding serine 33
Target	RPA32, Phospho (S33)
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P15927
NCBI	NP_002937.1
Dilution Range	WB - 1:10,000 - 1:25,000. Numerous bands representing proteins that migrate more slowly than 60 to 70 kDa have been observed during WB. Under some conditions, at least two of these may be more intense than the band that represents phospho RPA32-S33.
Expiration Date	12 months from date of receipt.



Detection of
human
Phospho
RPA32 (S33)
b...