



Product Datasheet

BRD8 Antibody (orb1530521)


Description	BRD8 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 225 and 275 of isoform 2
Target	BRD8
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9H0E9
NCBI	NP_631938.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 8 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human BRD8 by western blot....



Detection of human BRD8 by immunohistoch...



Detection of human BRD8 by western blot ...