


Product Datasheet

BRD2 Antibody (orb1530546)

Description	BRD2 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 550 and 600
Target	BRD2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P25440
NCBI	NP_005095.1
Dilution Range	WB - 1:1,000 - 1:10,000; IP - Not recommended
Expiration Date	12 months from date of receipt.



Detection of human BRD2 by western blot....