



Product Datasheet

H2AX Antibody (orb1530712)


Description	H2AX Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminus
Target	H2AX
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P16104
NCBI	NP_002096.1
Dilution Range	WB - 1:2,000 - 1:10,000; IHC - 1:40 to 1:200. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of mouse H2AX by western blot....



Detection of human H2AX by western blot....



Detection of human H2AX by immunohistoch...