



Product Datasheet

SMC2 Antibody (orb1530748)

Description	SMC2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 950 and 1000
Target	SMC2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O95347
NCBI	NP_006435.1
Dilution Range	WB - 1:1,000 - 1:10,000
Expiration Date	12 months from date of receipt.



Detection of human SMC2 by western blot....



Detection of human SMC2 by western blot....