


## Product Datasheet

### Myc Tag Antibody (orb1530881)

|                            |                                                                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Myc Tag Antibody                                                                                                                    |
| <b>Species/Host</b>        | Gallus                                                                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                                                                        |
| <b>Tested Applications</b> | ELISA, ICC, IP, WB                                                                                                                  |
| <b>Target</b>              | Myc Tag                                                                                                                             |
| <b>Preservatives</b>       | Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide                                                                       |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                           |
| <b>Concentration</b>       | 1 mg/ml                                                                                                                             |
| <b>Storage</b>             | 2 - 8°C                                                                                                                             |
| <b>Note</b>                | For research use only                                                                                                               |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                          |
| <b>Uniprot ID</b>          | <b>P01106</b>                                                                                                                       |
| <b>Dilution Range</b>      | WB - 1:1,000 - 1:30,000; IP - 1 - 4 µg/mg lysate; ICC - 1:100 - 1:400; ELISA - 1:1,000 - 1:30,000; for coating plates 1:100 - 1:500 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                     |



Detection of  
c-myc-  
tagged  
protein by  
wes...