


Product Datasheet

SARS-CoV-2 Spike RBD Antibody (orb1530887)

Description	SARS-CoV-2 Spike RBD Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	HEK293-derived SARS CoV-2 S-RBD protein Ala319-Phe541, with a C-terminal HIS6 tag
Target	SARS-CoV-2 Spike RBD
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P0DTC2
NCBI	YP_009724390.1
Dilution Range	WB - 1:10,000 - 1:25,000
Expiration Date	12 months from date of receipt.



Detection of coronavirus spike protein R...