

Product Datasheet

Sheep Albumin Antibody (orb1530999)

Description	Sheep Albumin Antibody
Species/Host	Rabbit
Reactivity	Sheep
Conjugation	Unconjugated
Tested Applications	Double Diffusion, IEP
Target	Sheep Albumin
Preservatives	Serum containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P14639
Dilution Range	IEP - Neat; DD - Neat
Expiration Date	12 months from date of receipt.