

Product Datasheet

Rabbit IgG-Fc Fragment Antibody (orb1531081)

Description	Rabbit IgG-Fc Fragment Antibody
Species/Host	Goat
Reactivity	Rabbit
Conjugation	AP
Tested Applications	ELISA, ICC, IHC, WB
Target	Rabbit IgG-Fc Fragment
Preservatives	50 mM HEPES pH 7.1, 0.1 M NaCl, 1 mM MgCl ₂ , 0.1 mM ZnCl ₂ containing 0.2% rAlbumin and 0.09% NaN ₃
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB - 1:5,000 - 1:50,000; IHC - 1:100 - 1:1,000; ICC - 1:100 - 1:1,000; ELISA - 1:5,000 - 1:50,000
Expiration Date	12 months from date of receipt.