

Product Datasheet

Rat IgG2b-Heavy Chain Antibody (orb1531215)

Description	Rat IgG2b-Heavy Chain Antibody
Species/Host	Goat
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Target	Rat IgG2b-Heavy Chain
Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	WB - 1:1,000 - 1:30,000; IHC - 1:200 - 1:2,000; ICC - 1:200 - 1:2,000; ELISA - 1:1,000 - 1:30,000; for coating plates 1:100 - 1:500
Expiration Date	12 months from date of receipt.