

Product Datasheet

GFAP Antibody (orb1532951)

Description	GFAP Antibody
Species/Host	Mouse
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant human alpha-helical GFAP fragment expressed in and purified from E. coli
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	P14136
Entrez	2670
Dilution Range	Immunofluorescence: 1:500; Immunohistochemistry: 1:500; Western Blot: 1:1,000
Expiration Date	12 months from date of receipt.