

Product Datasheet

CTNND1 / p120 Catenin Antibody (orb1533148)

Description	CTNND1 / p120 Catenin Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	LP-11111, synthetic peptide
Target	CTNND1 / p120 Catenin
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P (5 µg/ml)
Expiration Date	12 months from date of receipt.