

Product Datasheet

GLP1R / GLP-1 Receptor Antibody (N-Terminus) (orb1533271)

Description	GLP1R / GLP-1 Receptor Antibody (N-Terminus)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, IHC-Fr, IHC-P
Immunogen	Synthetic 16 amino acid peptide from N-Terminus extracellular domain of human GLP1R. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey (100%); Marmoset, Mouse, Rat, Sheep, Bovine, Hamster, Elephant, Rabbit, Pig (94%).
Target	GLP1R / GLP-1 Receptor
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IF, IHC, IHC-Fr, IHC-P (3 - 4 µg/ml)
Expiration Date	12 months from date of receipt.