

Product Datasheet

GLP1R / GLP-1 Receptor Antibody (Cytoplasmic Domain) (orb1533273)

Description

GLP1R / GLP-1 Receptor Antibody

Species/Host	Rabbit
Reactivity	Bat, Canine, Equine, Hamster, Human, Monkey, Mouse, Rabbit
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 16 amino acid peptide from 3rd cytoplasmic domain of human GLP1R. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset, Mouse, Dog, Bat, Hamster, Horse, Rabbit, Pig (100%); Sheep, Bovine, Elephant, Panda (94%); Rat, Opossum (88%); Turkey, Chicken (81%).
Target	GLP1R / GLP-1 Receptor
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 - 10 µg/ml)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558