

Product Datasheet

CHRM2 / M2 Antibody (Cytoplasmic Domain) (orb1533450)

Description

CHRM2 / M2 Antibody (Cytoplasmic

Species/Host

Rabbit

Reactivity

Bovine, Canine, Equine, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep

Conjugation

Unconjugated

Tested

IHC, IHC-P

Applications
Immunogen

Synthetic 16 amino acid peptide from 3rd cytoplasmic domain of human CHRM2. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Orangutan, Gibbon, Monkey, Marmoset, Mouse, Rat, Sheep, Dog, Bovine, Hamster, Elephant, Panda, Horse, Rabbit, Pig, Guinea pig (100%); Bat, Opossum, Turkey, Chicken (94%); Xenopus (88%); Platypus, Zebrafish, Stickleback (81%).

Target

CHRM2 / M2

Preservatives

PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide

Concentration

1 mg/ml

Storage

Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Dilution Range

IHC, IHC-P (6 - 15 µg/ml)

Expiration Date

12 months from date of receipt.