

Product Datasheet

SARS-CoV-2 Spike Glycoprotein Antibody (orb1533468)

Description	SARS-CoV-2 Spike Glycoprotein Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	LP-13540, synthetic peptide
Target	SARS-CoV-2 Spike Glycoprotein
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB
Expiration Date	12 months from date of receipt.