

Product Datasheet

SARS-CoV-2 Spike Glycoprotein Antibody (orb1533469)

Description	SARS-CoV-2 Spike Glycoprotein Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	LP-13541, synthetic peptide
Target	SARS-CoV-2 Spike Glycoprotein
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ELISA (10 µg/ml Optimal working dilution to be determined by the researcher), WB (Optimal working dilution to be determined by the researcher)
Expiration Date	12 months from date of receipt.