

Product Datasheet

ACE2 / ACE-2 Antibody (orb1533473)

Description	ACE2 / ACE-2 Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	LP-13632, synthetic peptide
Target	ACE2 / ACE-2
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ELISA (1:100), IHC-P (10 µg/ml), WB (1:1000)
Expiration Date	12 months from date of receipt.