

Product Datasheet

Rab5 antibody (orb153348)

Description

Goat polyclonal antibody to mouse Rab5. Rab5



Species/Host

Ra...

Goat

Reactivity

Bovine, Canine, Donkey, Equine, Feline, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep



Confocal
immunofluorescence
analysis of ...

Conjugation


Unconjugated

Tested
Applications

IF, IHC-Fr, IHC-P, WB

Immunogen

Purified recombinant peptide derived from within residues 115 aa to the C-terminus of mouse Rab5a, Rab5b and rab5c produced in E. coli.


Target

Rab5

Preservatives

PBS, 20% glycerol and 0.05% sodium azide

Form/Appearance

Polyclonal antibody supplied as a 200 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.

Western blot
analysis of spleen
cell lin...

Concentration

3 mg/ml

Storage

Store at -20°C for long-term storage. Store at 2-8°C for up to one month.

Note

For research use only

Application notes

The antibody solution should be gently mixed before use.

Isotype

IgG

Clonality

Polyclonal

Dilution Range

WB:1:250-1:2,000, IF:1:50-1:250, IHC-P:1:200-1:1,000, IHC-F:1:200-1:1,000

Expiration Date

12 months from date of receipt.