


Product Datasheet

Rab9 antibody (orb153350)

Description	Goat polyclonal antibody to mouse Rab9. Rab9 belongs to the small GTPase superfamily, Rab family. Th...
Species/Host	Goat
Reactivity	Bovine, Canine, Donkey, Equine, Feline, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant peptide derived from within residues 115 aa to the C-terminus of mouse Rab9a and Rab9b produced in E. coli.
Target	Rab9
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 200 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:250-1:2,000
Expiration Date	12 months from date of receipt.



Western blot analysis of Hepa cell line