

Product Datasheet

MAS1 / MAS Antibody (Extracellular Domain) (orb1533531)

Description	MAS1 / MAS Antibody (Extracellular Domain)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 17 amino acid peptide from 2nd extracellular domain of human MAS1. Percent identity with other species by BLAST analysis: Human (100%); Gorilla, Monkey (94%); Mouse (82%).
Target	MAS1 / MAS
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (20 µg/ml)
Expiration Date	12 months from date of receipt.