

Product Datasheet

ACAN / Aggrecan Antibody (Internal) (orb1533551)

Description	ACAN / Aggrecan Antibody (Internal)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic 19 amino acid peptide from internal region of human ACAN / Aggrecan. Percent identity with other species by BLAST analysis: Human, Monkey, Mouse, Rat (100%); Gorilla, Marmoset, Bat, Bovine, Dog, Hamster, Elephant, Panda, Horse (95%); Pig (89%).
Target	ACAN / Aggrecan
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P (2.5 µg/ml)
Expiration Date	12 months from date of receipt.