

Product Datasheet

MRGPRF Antibody (N-Terminus) (orb1533675)

Description	MRGPRF Antibody (N-Terminus)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 20 amino acid peptide from N-Terminus extracellular domain of human MRGPRF. Percent identity with other species by BLAST analysis: Human, Gorilla (95%); Chimpanzee, Monkey, Galago, Marmoset, Dog, Bovine, Elephant, Horse, Pig (89%); Orangutan, Mouse, Rat, Panda (84%).
Target	MRGPRF
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (11 - 17 µg/ml)
Expiration Date	12 months from date of receipt.