

Product Datasheet

T1R1 / TAS1R1 Antibody (N-Terminus) (orb1533681)

Description	T1R1 / TAS1R1 Antibody (N-Terminus)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IHC-P
Immunogen	Synthetic 18 amino acid peptide from N-Terminus extracellular domain of human TAS1R1. Percent identity with other species by BLAST analysis: Human (100%); Gorilla, Gibbon, Monkey (94%); Marmoset (83%).
Target	T1R1 / TAS1R1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	Flo, IHC, IHC-P (5 µg/ml)
Expiration Date	12 months from date of receipt.