

## Product Datasheet

### RORB / ROR Beta Antibody (Ligand-binding Domain) (orb1533771)

**Description**

RORB / ROR Beta Antibody (Ligand-

|                            |                                                                                                                                                                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                                                                                                       |
| <b>Reactivity</b>          | Bat, Bovine, Canine, Equine, Gallus, Hamster, Human, Monkey, Mouse, Rabbit, Rat, Sheep                                                                                                                                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                                                                                 |
| <b>Tested Applications</b> | IHC, IHC-P                                                                                                                                                                                                                                                                                                                   |
| <b>Immunogen</b>           | Synthetic 20 amino acid peptide from ligand-binding domain of human ROR Beta. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset, Mouse, Rat, Sheep, Bovine, Dog, Bat, Hamster, Elephant, Panda, Horse, Rabbit, Opossum, Chicken, Platypus (100%); Lizard (95%); Xenopus (85%). |
| <b>Target</b>              | RORB / ROR Beta                                                                                                                                                                                                                                                                                                              |
| <b>Preservatives</b>       | PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide                                                                                                                                                                                                                                                                               |
| <b>Concentration</b>       | 1 mg/ml                                                                                                                                                                                                                                                                                                                      |
| <b>Storage</b>             | Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.                                                                                                                                                                                  |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                                   |
| <b>Dilution Range</b>      | IHC, IHC-P (5 µg/ml)                                                                                                                                                                                                                                                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                                                                              |

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive&lt;br&gt;Research Triangle Park&lt;br&gt;Durham, North Carolina&lt;br&gt;27709. United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558