

Product Datasheet

RORB / ROR Beta Antibody (Hinge Domain) (orb1533787)

Description	RORB / ROR Beta Antibody (Hinge Domain)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P
Immunogen	Synthetic 15 amino acid peptide from hinge domain of human ROR Beta. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Marmoset, Sheep, Bovine, Dog, Bat, Panda, Horse, Rabbit, Opossum, Platypus (93%); Rat, Hamster, Elephant (87%); Mouse, Pig, Chicken, Lizard (80%).
Target	RORB / ROR Beta
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ELISA, IHC, IHC-P (20 µg/ml)
Expiration Date	12 months from date of receipt.