



Product Datasheet

LMO2 Antibody (orb153381)


Description	Mouse monoclonal antibody to LMO2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human LMO2 (AA: 1-158) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	18.4kDa
Entrez	4005
Dilution Range	RCM: 1/200 - 1/400; ELISA: Propose dilution 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Western blot analysis of human LMO2 usin...



Western blot analysis of HEK293 (Lane 1)...