



Product Datasheet

TWIST1 Antibody (orb153387)

Description	Mouse monoclonal antibody to TWIST1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human TWIST1 (AA: 9-74) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	21kDa
Entrez	7291
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Western blot analysis of human TWIST1 us...