



Product Datasheet

Mouse Lplunc1 (orb153389)


Description	Mouse monoclonal antibody to Mouse Lplunc1
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant fragment of mouse Lplunc1 (AA: 248-475) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	52.4kDa
Entrez	228801
Dilution Range	WB: 1/500 - 1/2000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Western blot analysis of mouse Lplunc1 u...



Western blot analysis of NIH3T3 (Lane 1)...